ApneaLink Air™ Tip Sheet



Below are some tips for your ApneaLink Air™ home sleep apnea test device.

- Call the Sleep Center at (313) 916-3222 if you have questions about the test or any
 equipment issues.
- Call the HFHS operator at (313) 916-3232 if you cannot reach the sleep center. The doctor on call will answer questions you have.
- Call 911 if you have a medical emergency during the test.

Tips to Use the Device:

- Review the ResMed device instructions that you got.
- Watch the YouTube video "How to use the ApneaLink Air™ Home Sleep Testing Device" by ResMed. You can get to the video by:
 - Search for the video on YouTube.
 - Type in this link: <u>youtube.com/watch?v=A1ujf-MyAsU</u>.
 - Scan the QR code.



To scan the QR code:

- 1. Open the camera on your smartphone.
- 2. Place the QR directly in the center of the camera and the link should appear on the screen.
- 3. This link will take you to the video.
- Keep anything that uses Bluetooth at least 15 feet away from your bed when you use this device. They can interfere with the sleep study. Some examples are:
 - Smartphones
 - Smart watches
 - Fitbits or other fitness trackers
 - Computers
 - Headphones
- Trim the nail that will have the finger probe on it. Remove any nail polish from that finger.
 - It's best to use a finger on your non-dominant hand.
- Use the tape you were given to secure the nose tubing and finger probe.
 - Tape the tube to your cheeks and close to your nose.
 - Wrap the tape around your finger where the line and finger probe connect. Tape the line to the top of your hand as well.
- You can go to the bathroom while you use the device. Do not take the device off or stop the study.
- Sleep as long as you can. We would like at least 7 hours of sleep data.

Pick-Up and Return the Device:

- 1. Go to the facility you are scheduled to pick-up or drop off your device.
- 2. Walk in to the sleep clinic desk. Let the Customer Service Representative know you are there to pick up or drop off a sleep device.

Columbus Sleep Center 39450 W. 12 Mile Rd. Novi, MI 48377 (248) 344-2460 3rd Floor Sleep Desk

Cottage Sleep Center 159 Kercheval Grosse Pointe Farms, MI 48236 (313) 640-2371 1st Floor Sleep Desk Garfield Road Sleep Center 43650 Garfield Rd. Clinton Township, MI 48038 (586) 203-1030 1st Floor Sleep Desk

Once you pick up your sleep device, you will have a questionnaire to fill out in MyChart. Please complete this as soon as you are able.

